

BISCAYNE ORTHODONTIC

LINE



The Biscayne Orthodontic Line provides simplistic functionality and a modern design, enabling the operator to work easily and efficiently while maintaining optimal patient comfort. With several accessory options and configurations available, it's easy to see how the Biscayne Orthodontic Line can enhance any practice.



With a power base and back, the Biscayne and Biscayne E.L. chairs provide the versatility to be used for either orthodontics or general dentistry. At a price point comparable to most fixed base chairs in the industry, the Biscayne and Biscayne E.L. are viable options for any practice.

CREATE YOUR OWN PACKAGE

*Possible configuration shown above

INSIDE SUMMIT DENTAL SYSTEMS

Since its founding in 1986, Summit Dental Systems (SDS) has been an internationally recognized manufacturer of high quality, durable dental equipment. Made in Pompano Beach, Florida, the SDS line of chairs provides the user with proven functionality, dependable mechanics, and undeniable value.

Constructed with steel, the Biscayne Chairs are designed to last a lifetime. With a quiet and smooth operation the chairs are comfortable, easy to use, and virtually maintenance free. Providing innovative solutions to any dental practice is our goal. We aim to sustain consistent quality in our dental equipment, and we stand by our mission by providing an industry leading five (5) year warranty.

The satisfaction of our customers is paramount, which is why our dedicated customer service team and highly skilled technical service department provide our customers with the support they need from their dental equipment manufacturer.



CREATE YOUR OWN PACKAGE

*Possible configuration shown above

CHOOSE YOUR ORTHODONTIC CHAIR

BISCAYNE

Chair Features:

- Hydraulic base and back
- Adjustable backrest speed
- Adjustable base down speed
- Automatic return
- Foot control
- Double articulating headrest
- Steel base and backrest assembly
- Powder-coat baked finish



BISCAYNE E.L.

Chair Features:

- Electro-mechanic base and back
- 1 Programmable pre-position
- Automatic return
- Foot control
- Double articulating headrest
- Armrests
- Steel base and backrest assembly
- Powder-coat baked finish
- Note: Recommended to be bolted to the floor



BISCAYNE FIXED BASE

Chair Features:

- Electro-mechanic back
- Foot control
- Double articulating headrest
- Steel base and backrest assembly
- Powder-coat baked finish
- Note: Recommended to be bolted to the floor



CHAIR OPTIONS



BISCAYNE SWIVEL ARMRESTS



BISCAYNE FIXED BASE ARMRESTS WITH ONE OUTWARD ROTATING ARM



MAGNETIC HEADREST



PONY TAIL HEADREST



V-SHAPED HEADREST

CHOOSE YOUR ORTHODONTIC CABINET



900SH

- Mobile cabinet
- 2HP automatic
- 3-way syringe
- SE & HVE
- Removable flat surface tray
- Water bottle system (option)
- Dimensions of cabinet top: 25³/₈" W, 18¹/₂" D



900SH WITH COMPUTER PACKAGE

- Mobile cabinet
- 2HP automatic
- 3-way syringe
- SE & HVE
- Monitor mount
- Keyboard drawer with mini keyboard
- CPU bracket
- Dual power outlet



900CWS WITH COMPUTER PACKAGE

- Mobile cabinet
- 2HP automatic
- 3-way syringe
- SE & HVE
- Monitor mount
- Keyboard drawer with mini keyboard
- Water bottle system
- Vented back door for computer tower storage
- Dual power outlet
- Dimensions of cabinet top: 26³/₈" W, 29¹/₂" D

Standard drawer paint color options:
white or gray

CABINET OPTIONS

- Fixed base
- Custom laminate
- Solid surface top
- Installation of custom monitor mount available


* Monitor not included with cabinets

BISCAYNE ORTHODONTIC



Standard Vinyl Color Chart

The colors and products in this brochure may vary from the actual materials. Ultra-leather and custom vinyl upholstery is available for an additional charge. Please see your authorized SDS dealer for accurate color samples and the most current product information.



Taupe



Bone



Cobalt



Wedgewood



Dusty Jade



Feather



Black



Desert Rose



800.275.3368 | Tel: 954.730.3636 | Fax: 954.730.3602

1280 SW 27th Avenue | Pompano Beach, FL 33069

www.summitdental.com