

PROUDLY PRESENTS THE

ViuluxLED II



Features & Benefits

- Cool white light increases color visualization
- Allows for more accurate diagnostics and color matching
- Light activation via motion sensor avoids cross contamination
- Precise lighting pattern
- Eliminates discomfort to patient's eyes
- Controllable brightness adjustment for more precision

Cutting edge lighting technology:

Light intensity: Variable 10,000 – 30,000 lux

LED life span: 50,000 hrs (25 yrs: 8 hrs/day)

Color temperature: 5,000 Kelvin

800.275.3368 | Tel: 954.730.3636 |

Fax: 954.730.3602 1280 S.W. 27th Ave |

Pompano Beach, FL 33309

www.summitdental.com

We proudly offer our products through our trusted authorized dealers