

6000EL Operatory Chair



The 6000EL Chair delivers quality and value ideal for hygienists and dentists. With a clean look and simple design, the 6000EL is practical for any dental setting. Economically priced without sacrificing quality, the 6000EL is a viable option for dentists who want a functional chair they can rely on for their practice. The 6000EL chair is standard without arms and can have arms added to the chair as shown below. As with all SDS chairs, the electro-mechanic chair is covered under our industry leading 5-year warranty.



SPECIFICATIONS

LIFT CAPACITY	440 lb (200 kg)
NET WEIGHT	230 lb (104 kg)
CHAIR WIDTH	26.5" (673mm)
CHAIR LENGTH WITH BACKREST DOWN	77" (1956mm)
MAX HEIGHT WITH BACKREST DOWN	34.5" (876mm)
MIN HEIGHT WITH BACKREST DOWN	20" (508mm)
MAX HEIGHT WITH BACKREST UP	57" (1448mm)
MIN HEIGHT WITH BACKREST DOWN	42" (1067mm)

FEATURES

- Steel base & backrest assembly
- Powder-coated baked finish
- Electro-mechanic base & back
- 1 programmable pre-position
- Automatic return
- Foot control
- Double articulating headrest
- Armrests Optional (one outward rotating arm)

6000EL

Operator and Hygiene Series

- A. 6000EL Chair
- B. Optional Armrests
- C. 1335 Light with Post
- D. 1402 2HP Automatic Delivery Unit
- E. 1615E Assistant's Arm with SE & HVE

The 6000EL Package is the perfect fit for all your dental or hygiene needs. Streamlined for success, the 6000EL Package can be equipped with the 1335 Light or Viulux LED II Light, the 2 or 3 handpiece delivery unit, and either the 1615E Assistant's Arm or the 1607 Cuspidor (both with SE and HVE). With these standard features and additional optional features, you'll see why the 6000EL Package is the economical solution to bring quality and functionality to your practice.



800.275.3368 | Tel: 954.730.3636 |
Fax: 954.730.3602
1280 SW 27th Ave | Pompano Beach, FL 33069
www.summitdental.com

