

# SDS Delivery Units

Over the Patient/Chair Mounted Options

## RELIABLE DELIVERY

Simplicity in design and maintenance has always been the SDS vision. The designs of our delivery units have been streamlined to be accessible and worry free. All our units are free of redundancies and complex features that can be costly and interruptive to a practice's day to day. The SDS units also provide the user peace of mind with durable ABS housing, steel frame, and a lifetime warranty on the control block diaphragm.



### 1440 EURO 4HP DELIVERY UNIT



\* SDS Chair Touchpad Available

#### FEATURES:

- 4 Handpieces Automatic
- Whip Arm Hosing
- Chair Mounted – Post or Swing Available
- 3-Way Air/Water Syringe
- Adjustable Water Coolant Controls
- Speed Controls for all Handpieces



### 1407D 4HP DELIVERY UNIT

#### FEATURES:

- 4 Handpieces Automatic
- Chair Mounted – Post or Swing Available
- 3-Way Air/Water Syringe
- Adjustable Water Coolant Controls
- Speed Controls for all Handpieces
- SDS Chair Touchpad Included



### 1407 3HP DELIVERY UNIT



\* SDS Chair Touchpad Available

#### FEATURES:

- 3 Handpieces Automatic
- Chair Mounted – Post or Swing Available
- 3-Way Air/Water Syringe
- Adjustable Water Coolant Controls
- Speed Controls for all Handpieces

### 1402 2HP DELIVERY UNIT

#### FEATURES:

- 2 Handpieces
- Chair Mounted – Post or Swing Available
- 3-Way Air/Water Syringe



800.275.3368 | Tel: 954.730.3636 | Fax 954.730.3602  
1280 SW 27<sup>th</sup> Ave | Pompano Beach, FL 33069  
[www.summitdental.com](http://www.summitdental.com)





# SDS Delivery Units

Over the Patient/Chair Mounted Options

## All-In-One Delivery Units

With a functional design and easy-to-use interface, the 1407H and 1402H delivery unit bring practicality to the operatory allowing for easy access to the handpieces and vacuum all located on the unit head. These delivery units are the perfect fit for any operator looking for a simple, reliable all-in-one 2HP unit for hygiene or general dentistry use.



### 1407H 2HP ALL-IN-ONE HYGIENE DELIVERY



\* SDS Chair Touchpad Available

#### FEATURES:

- 2 Handpieces Automatic
- Chair Mounted – Post or Swing Available
- SE, HVE, 3-Way Air/Water Syringe
- Adjustable Water Coolant Controls
- Speed Controls for all Handpieces



### 1402H 2HP ALL-IN-ONE HYGIENE DELIVERY UNIT



#### FEATURES:

- 2 Handpieces Automatic
- Chair Mounted – Post or Swing Available
- SE, HVE, 3-Way Air/Water Syringe
- Adjustable Water Coolant Controls
- Speed Controls for all Handpieces



800.275.3368 | Tel: 954.730.3636 | Fax 954.730.3602  
1280 SW 27<sup>th</sup> Ave | Pompano Beach, FL 33069  
[www.summitdental.com](http://www.summitdental.com)

