


## 9000PB/Palm Beach Chair Dimensions

(Includes Swing Arms Dimensions)


### With Base Fully Lowered

Note: Toe board moves approximately 8 inches forward when base is fully raised.



Distance from chair base with base fully lowered:

- A. Swing Light Arm 26"
- B. Swing Unit Arm 21"




Distance from toe board with base fully lowered:

- A. Swing Light Arm 17"
- B. Swing Unit Arm 12"


## With Base Fully Raised

Note: Toe board moves approximately 8 inches forward when base is fully raised.


Max Load Capacity  
**600 Lbs | 272.2 Kg**



### Footprint




## 9000PB/Palm Beach Junction Box Dimensions



## **Weight & Dimensions**

### **Shipping Information**



*Palm Beach Chair*  
477 Lb | 216 Kg